
Software Release Notes

CEMS Portal Version 5.0.1

This document provides the changes to CEMS version 5.01.

Purpose of the Release

The purpose of this document is to communicate major new features, changes and solved bugs in this release of CEMS software. It may also document workarounds and any feature restrictions.

Software Updates

- None

New Software Features

- None

Solved Software Bugs

- Clear alarm time format is not appropriate
- Hyperlink shows incorrect information for alarm/action level emails
- Halo dot and trend graph shows blank when data sent in decimal point
- Calibration SMS is not received by end users